Functional Plant Biology, 2014, 41, 223–226 http://dx.doi.org/10.1071/FP10900

PrometheusWiki Gold Leaf Protocol: gas exchange using LI-COR 6400

John R. Evans^{A,C} and Louis S. Santiago^B

^ABuilding 134, Research School of Biology, The Australian National University, Acton, ACT 0200, Australia. ^B3113 Batchelor Hall, University of California, Riverside, CA 92521, USA. ^CCorresponding author. Email: john.evans@anu.edu.au

Abstract. This leaf gas exchange protocol enables light or CO_2 response curves using a LI-COR LI-6400 portable photosynthesis system. This protocol originates in PrometheusWiki (http://prometheusWiki.publish.csiro.au/) where it has been tested and verified, and has received favourable user reviews. This reformatted and non-editable version is published as a Gold Leaf Protocol. For the most current version, including any user-commentary updates, readers may view the live version of the protocol at http://bit.ly/PWLicorGold.

Additional keywords: gas exchange, Li-Cor, photosynthesis, transpiration, leaf temperature, conductance.

Received 29 June 2010, accepted 28 October 2013, published online 11 February 2014

Introduction

This protocol explains how to make basic measurements with the LI-COR 6400 gas exchange instrument. The fluxes of CO_2 and water are measured on a portion of a leaf enclosed in a chamber where the CO_2 concentration, humidity, leaf temperature and light are controlled. The instrument calculates the rates of CO_2 assimilation and transpiration, stomatal conductance and the concentration of CO_2 in the intercellular airspaces within the leaf. Different leaves can be measured under a common ambient condition and responses to varying light or CO_2 can be determined to reveal information about the underlying biochemistry.

Units and abbreviations

Photo – photosynthetic rate, synonymous with CO₂ assimilation rate (μ mol CO₂ m⁻² s⁻¹); A – net CO₂ assimilation rate; Trmmol– transpiration rate (mmol H₂O m⁻² s⁻¹); Drierite – desiccant; A/C_i curve – CO₂ response curve; CO₂ concentration (μ mol CO₂ mol⁻¹ air); CO₂ partial pressure (μ bar or Pa); IRGA (infrared gas analyser); Irradiance (μ mol photons m⁻² s⁻¹); Reference air – air entering the leaf chamber; Sample air – air inside the leaf chamber; Cond – stomatal conductance to H₂O (mol H₂O m⁻² s⁻¹); VpdL leaf to air vapour pressure difference (kPa); PAR (photosynthetically active radiation, 400–700 nm); C_i – intercellular CO₂ concentration (μ mol CO₂ mol air⁻¹).

Materials and methods

Materials and equipment

- LI-COR 6400 with LED light source and CO₂ mixer (Note: This protocol is designed for versions 5 and 6(XT))
- CO₂ cartridges

- Drierite a desiccant
- Soda lime for CO₂ scrub
- Standard RS-232 cable to download data from LI-COR 6400
- LI-6400 software to download data
- LI-6400 Sim program to recompute data
- Tripod or field stand (optional, to mount clamp for reaching high leaves)
- PLUS permanent marker, razor blade, ziplock bag and cooler bag, for calculating leaf area
- LI-6020 battery charger

Operating the LI-6400

- (1) Put in batteries and attach sensor head to console (take care uncoiling the cable to prevent cable wind).
- (2) Unwind lid tension screw to tighten cuvette lid on sensor head.
- (3) Attach CO₂ cartridge: check O-ring is present and in good condition. Place cartridge where it lines up with O-ring, then screw on gently until it feels like it meets, then one last hard and fast screw turn to break the seal on the cartridge, make sure it is all the way on.
- (4) Turn on (using on/off switch).
- (5) In main menu, go to configuration menu, go to light source control, go to pick source and select the appropriate setting, press enter and escape to main menu.
- (6) Go to new measurements, line 2 and select temperature control. Set block temperature to control at your standard temperature.
- (7) In main menu, go to calibration menu, go to zero flow meter, you will have to wait ~10 s and then after verifying that the voltage is close to zero (otherwise manually adjust) you can hit OK (F5).

J. R. Evans and L. S. Santiago

- (8) (i) Humid conditions (where obtaining a valid H₂O zero may be difficult without dry N₂ gas). Make sure that there is fresh soda lime and drierite in the scrub columns. In the calibration menu, go to IRGA zeroing, make sure the CO₂ scrub is all the way on and that the desiccant is all the way off (bypass). Stay in this mode for 20 min to allow IRGAs to warm up. Record old zero values. Hit auto CO₂, but not auto H₂O or auto all. The H₂O should only be zeroed when running pure N₂ gas through the system. Verify good zero CO₂ values and record.
 - (ii) Dry conditions (replace used drierite in scrub column with fresh drierite). In the calibration menu, go to IRGA zeroing, turn the H₂O scrub on fully. Stay in this mode for 20 min to allow IRGAs to warm up. Verify good zero values and record. Hit auto H₂O and then turn the CO₂ scrub all the way on and the desiccant all the way off (bypass). Wait for a good zero to be reached (2–5 min) then verify good zero values and record. Hit auto CO₂, then exit calibration menu.
- (9) Go to new measurements menu. Disconnect leaf thermocouple and verify leaf temperature equals block temperature (row h). If not, adjust potentiometer underneath IRGA mounting bracket.
- (10) Go to menu (line) 1 and open the log file (F5), name the file (usually use date and species ID which is first 2 letters of each part of the scientific name, and any other information, like leaf # or plant #, e.g. 011607SOFA1–1) also put a short note, e.g. SOFA1–1
- (11) Go to menu (line) 2 and set the CO_2 mixer to reference 400 (menu > CO_2 mixer > reference >400) and the lamp to a reasonable value (1500 µmol m⁻² s⁻¹ for plants growing in sunlight, or 500 µmol m⁻² s⁻¹ for shade plants.) Decide on Temp (temperature) control if using batteries, this consumes power and requires more frequent battery changeover.
- (12) View the CO_2 trace (press]) and verify that CO_2 control has reached the set value. Match the IRGAs (line 1 F5). It is important to match the IRGAs every time CO_2 or H2OR is changed.
- (13) Check reference water (H2OR) is acceptable. If too low, add up to 10 mL water to the CO₂ scrub column (soda lime). Set flow rate if lower flow required.
- (14) Place the chosen leaf in the cuvette (Fig. 1), make sure to cover as much of the cuvette with the leaf as possible. Make sure the leaf is in contact with the thermocouple.
- (15) Go to menu (line) 3 and enter leaf area estimate if it does not fully cover the opening. Enter STOMRT (1 for equal density on top and bottom, 0 for stomata on one side only).
- (16) Once the cuvette is closed, watch the parameters (Photo, Cond, Ci, Tmmol). Ensure that transpiration does not lead to excessive humidity. Note: generally humidity greater than 85% is excessive. If condensation occurs, the gas exchange calculations will be wrong until the system is dried out.
- (17) View the graphs (press]) and when graph looks stable log the point using menu (line) 1, F1, or the button on the sensor head.
- (18) Light response curve or A-Ci curve. Go to menu 5, select AUTO PROG, pick type of curve, append to file (Y).



Fig. 1. Measuring leaf gas exchange in *Eucalyptus pauciflora* using a LI-6400. Photography by M. Harrison.

- (i) Light curve standard protocol: 1500, 1000, 500, 250, 120, 60, 40, 20, 10, 0), set min. wait time (usually 120 s), set max. wait time (usually 200 s), Y (yes) to match if change in CO_2 gets less than 20, Y (yes) stability definition is OK.
- (ii) CO₂ response curve: (*for* C_3 *species*) ambient, 50, 100, 150, 250, 350, 500, 700, 900, 1200, 700, 400. (*for* C_4 *species*) ambient, 30, 50, 80, 110, 150, 250, 350, 500, 700, 400. set min. wait time (usually 120 s), set max. wait time (usually 200 s), Y (yes) to match if change in CO₂ gets less than 200, Y (yes) stability definition is OK.
- (19) While the autoprograms are running, keep an eye on the following:
- display line C, Photo (photosynthesis) should be between 0–15 or higher, respiration in darkness 0 to –2, Cond (conductance) should always be positive, between 0.01–1, Ci should be positive (a negative value indicates a problem with leaf temperature or water calibration).
- display line D, VpdL, if greater than 2, then leaf is too hot, maybe too hot outside or too late in the day to continue (although OK for diurnals).
- display line G, can check PAR in (light source) and PAR out (light level outside cuvette)
- display line H, leaf temp should be slightly less than air temp
- display line I, battery charge, when below 10 V it will beep to change battery.
- (20) Can see graph of readings, menu 4, view graphs, (for XT software, press]), can scroll through but normally keep an eye on photosynthesis and conductance and check on stability, etc.
- (21) Can check on recorded points by going to menu 1, view file (F2), QuikPik Config (F1) (first time, afterwards can RePlot Graph), select light curve.
- (22) Also while light curve is running, there is time (it usually takes 20–30 min) to look for the next plant or do other measurements, such as leaf width, etc, just make sure

battery level is OK or that you can hear the warning beeps if it becomes low.

- (23) When curve is finished, the star by menu 2 will be gone (a star indicates an auto-program is running), then check the graph of recorded points (see step 21), if the graph is OK and you are done with the leaf, then menu 1, close file (f3).
- (24) Now go to sensor head, while holding the leaf in place, open the cuvette and mark the leaf coverage of the cuvette, using a thin sharpie permanent marker.
- (25) Cut leaf off cleanly at stem, using a razor blade and put in a ziplock bag, keep the bag cool and out of direct light (in a cooler) until you can put in a leaf area meter.
- (26) If another leaf is being measured, go back to 10, otherwise put the instrument in sleep mode, download data or turn off.

Download data from the LI-6400 (for XT machines, data can be removed using flashcard)

- (1) Back in the laboratory, you can download the data stored in the LI-6400, but first you must install the LI-6400 software.
- (2) Take the LI-6400 out of the case and put in the batteries, you do not need to attach sensor head, it will say 'Is IRGA attached?', select No (N), and enter, it will go to the normal main menu.
- (3) Connect the LI-6400 to the computer with a standard RS-232 cable.
- (4) Open the file exchange program on your computer.
- (5) On the LI-6400:
 - (i) from sleep mode, go Home in the utility menu (F5), select communications and then select file exchange mode
 - (ii) if LI-COR off, turn on, go to utility menu and select file exchange menu
- (6) On the computer hit the connect button to connect the computer to the LI-6400.
- (7) Drag and drop files from LI-6400 to the desired folder on the computer.
- (8) Once all files have been transferred, hit the disconnect button, exit computer program, turn off LI-6400 (same way as when you are taking measurements), unplug the cable.

Recompute data (after you know leaf area) from the LI-6400

- Open LI-6400 Sim program (this program acts as a simulator for the LI-6400, so open it and work in the same as you would on the real thing) (or recompute on the instrument before downloading). Note: version XT will provide results in an Excel file, meaning the LI-6400 Sim program is not essential. All formulae are already contained within the spreadsheet (for recomputing results based on leaf area etc). [See Appendix 1 for example data for CO₂ response curves and a light response curve.]
- (2) Open the utility menu, choose recomputed data, Yes (Y) open window

- (3) On the file page, open the specific file, choose the data #1 tab and the data will be displayed in table form
- (4) In the column titled area, enter the correct leaf area (if less than 6, if 6 or greater then no need to recompute) down the column, then hit the recompute button (a triangle button that looks like a play button)
- (5) Save this file as the same name but with RCMP (recomputed) added
- (6) Open the next file and repeat until finished
- (7) Close out of program

Turning off the LI-6400

- (1) Turn off at switch
- (2) Loosen the cuvette lid to prevent gasket compression during storage
- (3) Remove the batteries and recharge
- (4) Fill in log book
- (5) Put everything away in the case.

Notes and troubleshooting tips for LI-6400

Synonyms: LI-COR 6400, LI-6400, 6400; protocol steps will differ slightly depending on software version being used. *Battery*: When battery is low, the machine will beep slowly and then faster as you run out of time. Get fresh batteries, remove one old battery and replace it. Then do the other, being careful not to be too hasty or to hit the off switch, as this will stop the instrument and you will have to start all over again.

Links to resources and suppliers

LI-6400/LI-6400XT Instruction Manuals and Technical Notes LI-COR Biosciences sell a photosynthesis educational resource package that is a series of PowerPoint lectures covering both theory and practice for multiple measurement protocols. It includes guidelines on what environmental control settings to select and why as well as what is being measured and how. Part number is 6400–520. <u>LI-COR Biosciences – Photosynthesis</u> <u>Educational Resource Package</u>. A practical guide to using the LI-6400 is provided in the Supplementary Material, which is available online.

Useful references

- Farquhar GD, von Caemmerer S, Berry JA (1980) A biochemical model of photosynthetic CO₂ assimilation in leaves of C₃ species. *Planta* 149, 78–90. doi:10.1007/BF00386231
- Ögren E, Evans JR (1993) Photosynthetic light response curves. 1. The influence of CO₂ partial pressure and leaf inversion. *Planta* **189**, 182–190. doi:10.1007/BF00195075
- Sharkey TD, Bernacchi CJ, Farquhar GD, Singsaas EL (2007) Fitting photosynthetic carbon dioxide response curves for C-3 leaves. *Plant, Cell & Environment* **30**, 1035–1040. doi:10.1111/j.1365-3040.2007. 01710.x
- von Caemmerer S, Farquhar GD (1981) Some relationships between the biochemistry of photosynthesis and the gas exchange of leaves. *Planta* 153, 376–387. doi:10.1007/BF00384257

3

4

5

6

41.2

25.3

9.62

3.83

0.31

0.234

0.17

0.133

159

207

298

345

1001

499

199

100

204.357

121.648

52.069

26.655

(a) Wheat		Photo	Cond	Ci	Trmmol	VpdL	Tair	Tleaf	TBlk	CO2R	CO2S	H2OR	H2OS	PARi
	1	24.6	1.01	323	8.18	0.961	20.89	22.07	19.98	401.31	377.22	10.961	17.949	1500
	2	0.773	1.03	46.4	8.23	0.954	20.89	22.09	19.98	49.67	48.65	11.038	18.071	1499
	3	5.35	1.04	83.1	8.25	0.948	20.88	22.07	19.98	100.13	94.79	11.054	18.106	1500
	4	13.7	1.06	158	8.29	0.939	20.88	22.04	19.98	200.21	186.91	11.059	18.14	1500
	5	20	1.07	240	8.31	0.933	20.89	22.03	19.98	301.38	281.92	11.086	18.183	1499
	6	23.9	1.08	326	8.31	0.93	20.89	22.02	19.98	400.19	376.64	11.096	18.198	1500
	7	26.6	1.09	516	8.35	0.928	20.89	22.02	19.98	600.59	573.27	11.088	18.223	1500
	8	26.7	1.1	711	8.36	0.924	20.89	22.01	19.98	800.42	771.59	11.104	18.244	1499
	9	27.9	1.1	903	8.4	0.93	20.89	21.99	19.98	1000.78	969.44	10.97	18.149	1499
(b) Paspalu	m dilatati	um												
	10	4.7	0.396	16.9	4.28	1.09	25.37	25.77	25.04	41.17	37.96	20.803	23.566	1501
	11	5.19	0.399	24.2	4.41	1.12	25.36	25.9	25.04	50.96	47.4	20.759	23.594	1501
	12	11	0.387	41.6	4.38	1.14	25.36	26.02	25.04	99.72	92.2	20.738	23.559	1501
	13	18.8	0.392	102	4.38	1.13	25.36	25.93	25.04	200.46	187.52	20.699	23.519	1501
	14	22.8	0.38	274	4.28	1.14	25.36	25.9	25.04	400.71	384.59	20.647	23.407	1501
	15	24.8	0.295	429	3.57	1.2	25.35	25.97	25.04	600.35	582.58	20.581	22.884	1501
	16	25.3	0.234	586	2.99	1.25	25.36	26.02	25.04	800.34	782.09	20.552	22.483	1501
	17	26.3	0.154	681	2.16	1.34	25.37	26.21	25.04	1001.03	982.26	20.515	21.91	1501
	18	23.5	0.185	166	2.5	1.3	25.37	26.14	25.04	399.49	383.32	20.552	22.165	1501
(c) Maize	light res	ponse curv	ve with n	nultiple f	lashes for	the fluor	escence	saturatio	n pulse					
Obs	-	Photo		Cond	Ci		Fm'		Fs	PhiPS2		ETR		PARi
Multiple fla	sh													
	1	54.6		0.364	124		583.7		389.5	0.333		282.926		2001
	2	49.4		0.35	141		654.5		389.6	0.405		257.816		1499

756.9

904.9

982.7

1011.3

393.4

386

379

376.2

0.48

0.573

0.614

0.628

Appendix 1. Example data for CO₂ response curves for (a) wheat (C₃) and (b) paspalum (C₄) leaves and (c) a light response curve for maize (C₄)

